

BYZ-810 Double-channel Serial Syringe Pump



BYZ-810T

Injection modes: Simple mode.

BYZ-810TU

3 Injection modes are available: Simple mode, Vol/T mode, Vol/W mode;
Human voice alarm system;
History record: up to 2000pcs records;
Drug library.

Brand-new Appearance Better Sensibility

> Design as portable



- Appearance is upgraded to provide better sensibility;
- Double bigger display windows enable more legibility and convenience.

Intelligent Safety More Care

- Compatible with all brands of syringes
- VTBI (Volume to be injected) is available.



- Intelligently identify syringe size: 10ml, 20ml, 30ml, 50ml (60ml);
- Safer reassurance is achieved by alarm combination of human voice, legible words and lights;
- Comprehensive drug library is able to satisfy various clinical needs.



$\leq \pm 3\%$

Precision

$\frac{0.1\text{ml/h}}{1500\text{ml/h}}$

Rate

0.1~5.0ml/h

KVO RATE

≥ 1000

Drug library



Features



Dual CPU

Guarantee for safe work.



Dynamic Pressure System

(DPS) 3-level occlusion pressure, dynamic safety monitoring




System Self-test

A system self-test is conducted every time it is turned on, including checks of the alarm indicator, volume, etc., providing a safeguard for each use.



Double Display Windows



High Accuracy



Independent Control Center