



Talking Pet Infusion Pump



BEYOND Vet Infusion Pump

Model: BY-100I

- Unique human voice alarm system
- Easy to operate with infusion memory function
- Free setting of infusion volume
- Compatible with universal brands of Iv sets
- Multiple infusion modes
- Dual CPU to make sure infusion safety



Professional



Honest



Simple



Efficient

Technology Serves Human Health

Vet Infusion Pump

 Model: **BY-100I**

Infusion rate range	1ml/h ~ 1200ml/h (Increment 1ml/h)
	1d/min ~ 266d/min
Infusion accuracy	ml/hMode, within $\pm 5\%$; d/min, within $\pm 3\%$
Mechanical accuracy	within $\pm 2\%$
FF rate	100 ml/h ~ 1000 ml/h , Adjustable (Increment 100ml/h)
Infusion volume	1ml ~ 9999ml
Cumulative infusion Volume	0ml ~ 9999.9ml (Increment 1ml/h)
KVO rate	1ml/h ~ 10ml/h (Increment 1ml/h)
Alarm	Infusion occlusion,Door open, Infusion set mismatch, Battery empty,Battery low, Air bubble,forgotten operation, Dropsensor is pulled out,Near infusion finished, Etc.
Air bubble sensor	Ultrasonic detection method, with a detection sensitivity of $\geq 25\mu\text{L}$
Power source	Two options of power supply: 100~240VAC,50/60Hz; Built-in Li battery 11.1V, The rechargeable lithium-ion battery pack has a capacity of $\geq 1800\text{mAh}$. After a 12-hour charge with a new battery, it can support the infusion pump to operate continuously for more than 6 hours at a rate of 25ml/h
Iv set	Compatible with universal brands of Iv sets(including light-proof sets)

 **24^{hour} service hotline** // +86-731-82564299 / +86-19807315856 (WeChat & WhatsApp)

Email : sales@beyondmedtech.com URL: www.beyondmedtech.com

Add: Beyond Zone, Zhixian Rd, Xueshi Street, Yuelu District, 410208, Changsha, China

Version number: CY-02-47-02 A/0-002



Beyond Website