

BYOND

RMP Series Portable Outlet Power Supply



Reliable power for your ResPlus, ResFree, or ResHope NIV.
Move freely. Stay powered.

Technology Serves Human Health

Applications:



Patient Transport



Short Trips



Emergency Backup for Power Outages



Outdoor Camping or Vehicle Backup



Key Features

- Functions as an uninterruptible power supply (UPS) for home use
- Equipped with DC, USB, and Type-C ports to support a wide range of devices and charging needs.
- Provides pure sine wave power for safe use with ResPlus, ResFree, or ResHope CPAP, BiPAP, or ventilator machines.
- Silicone anti-slip pad for enhanced stability on various surfaces.
- Soft night light for added convenience during nighttime use.



Parameters

Model	RMP100	RMP200
Capacity	92.5Wh	249.7Wh
Run Time on Full Charge	Powered for 3.5 hours with heated tube (20°C) on standard mode.	Powered for 8 hours with heated tube (20°C) on standard mode.
Output	1) Supports a maximum unidirectional DC-24V 3A output (via the DC-6330 port on the back). 2) Supports Type-C/PD-18W unidirectional output (supports 5V/9V/12V output). 3) Supports USB-A/QC-18W unidirectional output (supports 5V/9V/12V output). When three ports are used simultaneously, the DC port must be prioritized to ensure a maximum output of 72W. The combined maximum output for the USB and Type-C ports will be 18W.	
Operating Altitude	70kPa to 106kPa $\approx$ 2824m (9265 ft.) above sea level.	
Dimensions	250*155*30mm	
Weight	815g	
Safety Protection	Comprehensive protection includes temperature protection and short-circuit protection. Over-current protection is rated at 15A-20A. Battery over-discharge protection voltage is 13.75V $\pm$ 0.5V.	

 **24^{hour} service hotline** // +86-731-82564299 / +86-19807315856 (WeChat & WhatsApp)

Email : sales@beyondmedtech.com URL: www.beyondmedtech.com

Add: Beyond Zone, Zhixian Rd, Xueshi Street, Yuelu District, 410208, Changsha, China

Version number: CY-02-43-03 A/0-002



Beyond Website