

**BYOND**



## Disposable Syringes for Pumps

Precise Delivery | Safe & Reliable | Full Range of Specifications



Scan for video

Technology Serves Human Health

## + Versatile Sizes

Available in sizes ranging from 2ml to 60ml, these syringes are designed to meet a wide spectrum of clinical requirements—from precise small-dose administration to high-volume infusions.



## + Light-Proof Option

Specially designed for photosensitive medications, the light-shielding material maintains drug stability in sensitive environments.



## + Smooth Operation

Optimized resistance and precise matching between barrel and plunger provide a seamless push experience, perfect fit for high-precision syringe pumps with stable, uninterrupted delivery.



## + Secure Threaded Connection

Screw-type interface ensures a firm lock with the pump connector, preventing detachment or loosening during infusion for maximum safety.



| Syringe | Model | Packaging     | Connector | Light-Protective |
|---------|-------|---------------|-----------|------------------|
|         | 2ml   | 1800 Pcs/Case | Slip/Lock | Optional         |
|         | 3ml   | 1800 Pcs/Case | Slip/Lock | Optional         |
|         | 5ml   | 1200 Pcs/Case | Slip/Lock | Optional         |
|         | 10ml  | 800 Pcs/Case  | Slip/Lock | Optional         |
|         | 20ml  | 600 Pcs/Case  | Slip/Lock | Optional         |
|         | 50ml  | 300 Pcs/Case  | Slip/Lock | Optional         |

