



EaseFit FM3 Face Mask

S

M

L

Technology Serves Human Health

Features

Light Weight Upgrade

The newly upgraded technology balances daily comfort with lightness. The lighter structure reduces facial pressure and meets the convenient needs of daily respiratory treatment. At the same time, it is durable and an economical choice for home treatment.



Safety and Quiet

Anti-clogging micro exhaust hole design makes it safer for users to wear. Scientifically distributed exhaust design greatly reduces the airflow noise to below 35 decibels, allowing you to enjoy high-quality sleep moments.




Soft Fit

New upgraded liquid silicone technology, adaptive silicone pads fits the facial curve, gently wraps, reduces the pressure on the bridge of the nose, and adapts to various face shapes.



Skin-Friendly and Breathable

The headgear uses highly breathable and skin-friendly materials to reduce heat and friction, improve user experience.



Easy Assembly

It is fully upgraded to a quick-release buckle to simplify the assembly process, making it easier to wear. Combined with lightweight materials, it improves user comfort and convenience.



24 hours service hotline // +86-731-82564299 / +86-19807315856 (WeChat & WhatsApp)

Email : sales@beyondmedtech.com URL: www.beyondmedtech.com
Add: Beyond Zone, Zhixian Rd, Xueshi Street, Yuelu District, 410208, Changsha, China
Version number: CY-02-57-01 A/0-002



Beyond Website