



EaseFit FMI Face Mask

S

M


L

Technology Serves Human Health

Features

Intelligent Adjustment

Independent invention patent of forehead stepless knob design can adjust the angle more flexibly and easily according to the different face shapes of users. It can achieve precise fit by simple operation.



Safety and Quiet

Anti-clogging micro exhaust hole design makes it safer for users to wear. Scientifically distributed exhaust design greatly reduces the airflow noise to below 35 decibels, allowing you to enjoy high-quality sleep moments.




Soft Fit

A double-layer air rebound silicone pad could reduce the pressure on the skin, gently wrap and make it soft and traceless after long-term wearing, suitable for a variety of face shapes.



Skin-Friendly and Breathable

The headgear use highly breathable and skin-friendly materials to reduce heat and friction, improve user experience.



Light and Convenient

The full face mask provides more comprehensive airflow support and is suitable for people with nasal obstruction, mouth breathing, severe snoring, and chronic obstructive pulmonary disease. Lightweight design reduces facial burden and is comfortable to wear for a long time. The rotating quick-release connector supports 3-second disassembly, simplifying daily wear and cleaning maintenance.



24 hours service hotline // +86-731-82564299 / +86-19807315856 (WeChat & WhatsApp)

Email : sales@beyondmedtech.com URL: www.beyondmedtech.com
Add: Beyond Zone, Zhixian Rd, Xueshi Street, Yuelu District, 410208, Changsha, China
Version number: CY-02-09-01 A/0-002



Beyond Website