



Veterinary Infusion Pump

New Design

- Exquisite and mini, portable
- 3.5 Inch large colorful touch-screen and key-board, dual control methods.



Simple and Convenient

- Auto-lock Screen**
The lock screen duration can be adjusted in multiple levels; prevention for operation errors.
- Night Mode**
Design for night time safe operation.
- History Record Cassette**
Up to 50000pcs records.
Trackable treatments are ensured.
- System Self-test**
System self-test is conducted every time turning on, safeguard for every-time use.

Accurate and Safe



- Large capacity lithium battery
- IP44 High-level dust-proof and waterproof
- Adjustable multi-level bubble detection sensitivity to minimize air embolism risk
- Multiple intelligent and intuitive alarms
- Nighttime tubing illumination function
- DPS intelligent occlusion management

Configuration

Parameter	Specification			
Model	VP6 Vet	VP6 Vet (Air)	VP6 VET (PRO)	VP6 VET (PLUS)
Applicable Infusion Sets	1ml/15d, 1ml/19d, 1ml/20d, 1ml/60d			
Flow Rate	ml/h Mode: 0.1 ~ 1200 ml/h (in 0.01ml/h increments) Drops/min Mode: 1 (400ml/h*drop/60) d/min, (in 1 d/min increments)			
WIFI	N/A	YES	N/A	YES
Heating Tube	N/A	N/A	YES	YES
Pressure Level	5 Level			
VTBI	0.1~9999.99ml (in 0.01ml increments)			
Cumulative Infusion Volume	0~9999.99ml (in 0.01ml increments)			
Rate (ml/h) Accuracy	ml/h mode:≤±5%(after calibration), Drop mode:≤±5%(after calibration)			
KVO Rate (manual mode)	0.1~30.0ml/h (in 0.01ml/h increments)			
Drug Library	Up to 6000 drugs can be stored			
Dimensions	132mm(L) x 101mm(W) x 164mm(H)			
Weight	About 1.4kg			