



Sunfusion® PI3

Central Infusion Station

- + Multiple Channels
- + Power Integration

Technology Serves Human Health

Sunfusion® PI3 System Composition



Top Cover



Plug-in box



Sunfusion® Semi-1



Sunfusion® Semi-2



Sunfusion® Erch-3
(WIFI connectivity)



Sunfusion® Erch-6
(WIFI connectivity)

Integration—Save the space

+ Power Integration

- Only one power cable is required for charging and powering all pumps in the infusion workstation.
- Super design for tube arrangement and fixing of drop sensor in the two sides of plug-in box.



+ Combined Flexibly

- Modular design, free combination, small size, light weight and easy transportation.
- Each slot with auxiliary pallet, insert the pump through the rail, easy to operation.
- Stackable and stacked track design, easy to stack pumps without any tools.

Multiple Channels—Safe, exible and convenient

+ Hot Plugging Design

- Innovation hot plugging design, can pull out the pump and keep infusion when patient leaving bed.
- Infusion module, high performance and safe.



+ Combined Stacking

- plug-in and modular design: 2-12 pumps can be combined flexibly as required.

24hour service hotline // +86-731-82564299 / +86-19807315856 (WeChat & WhatsApp)

Email : sales@beyondmedtech.com URL: www.beyondmedtech.com

Add: Beyond Zone, Zhixian Rd, Xueshi Street, Yuelu District, 410208, Changsha, China

Version number: CY-02-55-02 A/0-002



Beyond Website