



Sunfusion® WI3

Multi-channel Infusion Dock-Station

- + Unified Audio-Visual Alarms
- + Multiple Channels
- + Power Integration

Technology Serves Human Health

BEYOND

Sunfusion® WI3 System Composition



Alarm controller



Plug-in box



Sunfusion® Semi-1



Sunfusion® Erch-3
(WIFI connectivity)



Sunfusion® Semi-2



Sunfusion® Erch-6
(WIFI connectivity)

Integration—Save the space

+ Power Integration

- Only one power cable is required for charging and powering all pumps in the infusion workstation.
- Super design for tube arrangement and fixing of drop sensor in the two sides of plug-in box.



+ Combined Flexibly

- Modular design, free combination, small size, light weight and easy transportation.
- Each slot with auxiliary pallet, insert the pump through the rail, easy to operation.
- Stackable and stacked track design, easy to stack pumps without any tools.



Multiple channels—Safe, flexible and convenient

+ Hot Plugging Design

- Innovation hot plugging design, can pull out the pump and keep infusion when patient leaving bed.
- Infusion module, high performance and safe.



+ Combined Stacking

- plug-in and modular design: 2-12 pumps can be combined flexibly as required.



Audible and Visible Alarm

+ Alarm Function

- The multi-channel infusion dock-station can provide the alarm status of infusion pump / syringe pump to users. Once the system triggers an alarm or other abnormal conditions, the alarm controller will receive the news and generated Audible and Visible Alarm to ensure the safety of infusion.



24 hour service hotline // +86-731-82564299 / +86-19807315856 (WeChat & WhatsApp)

Email : sales@beyondmedtech.com URL: www.beyondmedtech.com

Add: Beyond Zone, Lijiacun Rd, Xueshi Street, Yuelu District, 410208, Changsha, China

Version number: CY-02-55-01 A/0-002



Scan for video