

BYOND



Sunfusion[®] Series

Enteral Feeding Pump

ENP-1

ENP-5

Technology Serves Human Health

Sunfusion[®] Series Enteral Feeding Pump

Precise Feeding and Simple Operation

+ Simple design, easy operation

- || 3.5-inch touch display integrated, information is clear at a glance.
- || Dual operation of knob and touch screen makes it more convenient to use.
- || Equipped with a detachable battery compartment for easy replacement.



+ Precise feeding, safe protection

- || Built-in high-performance imported motor, feeding error $\leq 5\%$.
- || Graphical voice alerts, intuitive and efficient for safety.
- || Multiple working modes: continuous, intermittent, and reverse mode.
- || DPS intelligent blockage management helps medical staff manage infusion line blockage problems.

Smart Heating and Information Management

+ Intelligent heating, worry-free care

- || Built-in heating system, dual CPU temperature control and monitoring, to achieve automatic constant temperature heating. (ENP-5 only)



+ Information management, efficient and convenient

- || Optional WIFI function allows easy data transmission to the infusion central station software, real-time monitoring of feeding and improved medical work efficiency.

$\leq \pm 5\%$

Flow Accuracy

$\frac{1\text{ml/h}}{1200\text{ml/h}}$

Feeding Speed

IP34

Waterproof Rating

1ml/h

Minimum Step

1-99999.9ml

Preset Volume Range

24^{hour} service hotline // +86-731-82564299 / +86-19807315856 (WeChat & WhatsApp)

Email : sales@beyondmedtech.com URL: www.beyondmedtech.com

Add: Beyond Zone, Zhixian Rd, Xueshi Street, Yuelu District, 410208, Changsha, China

Version number: CY-02-28-04 A/0-002



Beyond Website