

BYOND




Sunfusion[®] Series

Nutrition Pump

ENP-1

Technology Serves Human Health

Sunfusion[®] Series Nutrition Pump



+ Intelligent Automatic

- || Enables micro-dose, continuous, and timed delivery control.
- || Real-time cumulative feeding volume display supports clinical observation.
- || Optional Wi-Fi connectivity improves the efficiency of medical monitoring.
- || Auto language switching (Chinese/English) for seamless operation.

+ Safe and Worry-Free

- || Highly safe medication administration, with a lithium battery-assisted power supply ensuring continuous medication delivery even during AC power outages.
- || High-tech intelligent alarm system with animated instructions.
- || Multiple adjustable pressure levels and real-time dynamic pressure display.
- || Convenient data access with historical record storage.
- || Highly safe device with IP34 dust and water resistance.



+ User-Friendly Design

- || Compact, lightweight, and easy to carry.
- || Fully transparent pump door design allows for clear visualization of nutrient delivery.
- || Multiple operating modes: continuous, intermittent, and reverse pumping.
- || 3.5-inch high-definition display for quick access to key information.
- || Dual-control input via touchscreen and keyboard, plus a jog dial for effortless operation.
- || Removable battery compartment for easy replacement.

≤ ±5%

Flow Accuracy

1ml/h
1200ml/h

Feeding Speed

IP34

Waterproof Rating

1ml/h

Minimum Step

1-99999.9ml

Preset Volume Range

24 hours service hotline // +86-731-82564299 / +86-19807315856 (WeChat & WhatsApp)

Email : sales@beyondmedtech.com URL: www.beyondmedtech.com

Add: Beyond Zone, Zhixian Rd, Xueshi Street, Yuelu District, 410208, Changsha, China

Version number: CY-02-28-02 A/0-002



Beyond Website